

# MAT 100 MEMORY FOAM

## MAGNETOTHERAPY MAT

MEMORY FOAM MAT FOR TREATMENTS WITH PULSED ELECTROMAGNETIC FIELDS. IT CAN BE USED WITH ALL THE MAGNETOTHERAPY DEVICES OF THE GLOBUS MAGNUM LINE.

MAT 100 IS IDEAL FOR TREATMENTS ON LARGE SURFACES, LIKE THE VERTEBRAL COLUMN.

THE THREE-LAYER PADDING IS VERY COMFORTABLE FOR THE TREATED PART AND THIS WILL BE PARTICULARLY USEFUL IN CASE OF STRONG PAIN AND INFLAMMATION.

### TECHNICAL SPECIFICATIONS

#### Size

70x44x6,5 cm

#### Internal padding

Double memory-foam surface layer with different density and polyurethane layer to insert the solenoids

#### External coating

Washable cover in hypoallergenic fabric with slip-resistant underside

#### Internal coating

Washable cover in water-resistant, antibacterial fabric

#### Number of solenoids

4

### RECOMMENDED IN CASE OF


Back pain and dorsalgia

Arthrosis and arthritis

Muscle pain



### INTERNAL SECTION OF THE DIFFERENT LAYERS OF THE MAT



# MAT 100 MEMORY FOAM

MAGNETOTHERAPY MAT

Rev. En\_01\_2022