

MAGNUM 3000 PRO

PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

POWERFUL PROFESSIONAL MAGNETOTHERAPY EQUIPPED WITH MANY SPECIAL PROGRAMS AND FUNCTIONS THAT MAKE IT EFFECTIVE THE TREATMENT OF SEVERAL OSTEOARTICULAR PATHOLOGIES AND THE RICOVERY AFTER ORTHOPEDIC SURGERY. POSSIBILITY OF LOCAL AND TOTAL-BODY TREATMENTS.

TECHNICAL SPECIFICATIONS

Display

Size of the visible area 2,6"

Channels

2 independent

Intensity can be emitted with flexible diffuser

- 200 Gauss on average for channel
- 400 Gauss on average on the whole (for 2 channels)

Frequency

5-200Hz

Duty Cycle

5-50%

Free programming

Possibility of adding up to 10 custom programs

1+1 function

Possibility of running 2 different programs at the same time, using 1 channel each

Power supply


- Mains
- Rechargeable batteries (not included)

Size | Weight

170x220x60 mm | 1152 gr

STANDARD EQUIPMENT

- 1 Bag
- 1 MAGNUM 3000 PRO magnetotherapy device
- 2 Flexible solenoids 300x100 mm
- 1 Magnum Test
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- MAGNUM 3000 PRO with 2 flexible solenoids
- MAGNUM 3000 PRO with 2 soft solenoids
- MAGNUM 3000 PRO with 2 solenoids Pocket Pro



CE
0476



EASY PROGRAM



1+1 MODE



LAST 10



PROGRAMMABLE



RENTAL OPTION



RUN TIME / EDIT PARAMETERS



MAINS SUPPLY

GLOBUS
ITALIAN EXCELLENCE

MAGNUM 3000 PRO

PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

PROGRAM LIST

ARTHROSIS/ARTHRITIS

Knee arthrosis

Cervical arthrosis

Arthritis

Arthrosis

Hand arthrosis

Hip arthrosis

Shoulder arthrosis

Column arthrosis

PAIN

Algodystrophy - Sudek disease

Back pain

Joint pain

INFLAMMATIONS

Anti-inflammatory

Tendinitis

Epicondylitis

Epitrochleytis

Carpal tunnel

Shoulder tendinitis

ORTHOPEDICS

Post ACL surgery

Cartilage lesions

Knee posthesis

Hip prosthesis

Spinal fusion

Osteotomy

BONE PATHOLOGIES

Fractures

Scaphoide fracture

Wrist fracture

Shoulder fracture

Collarbone fracture

Forearm fracture

Pelvic fracture

Tibial fracture

Femoral fracture

Delayed fracture consolidation

Fracture with plaster

Osteonecrosis

Osteoporosis

Pseudoarthrosis

Bone edemas

TRAUMAS / EDEMAS / WOUNDS

Ankle sprain

Edemas- contusions

Ulcers

NIGHT PROGRAMS

Night program 15Hz

Night program 50Hz

Night program 75Hz

DEFAULT PROGRAMS WITH DIFFERENT FREQUENCY

Multi-frequency program

Program 5Hz

Program 10Hz

Program 15Hz

Program 20Hz

Program 25Hz

Program 30Hz

Program 35Hz

Program 40Hz

Program 45Hz

Program 50Hz

Program 55Hz

Program 60Hz

Program 65Hz

Program 70Hz

Program 75Hz

Program 80Hz

Program 85Hz

Program 90Hz

Program 95Hz

Program 100Hz

Program 120Hz

Program 140Hz

Program 160Hz

Program 180Hz

Program 200Hz